



AMSD Alignment Lessons Learned

Patrick J. Reardon, James B. Hadaway
The University of Alabama in Huntsville

H. Philip Stahl, Ron Eng, Jay Carpenter, &
Jeff Kegley
Marshall Space Flight Center

September 16, 2003



Outline



- Overview of Test Procedure
- Selected Details
- Repeatability and Reproducibility Tests
- Summary



Overview of Test Procedure



- Phase I: Ambient Pressure
 - Align Instantaneous Phase Interferometer (IPI) to Null
 - ◆ Record position of null in pixel space
 - Rough Alignment of surface to IPI/Null
 - ◆ Using spot projectors
 - ◆ Gets surface within 0.5-1.0mm of alignment
 - ◆ Requires fiducials to be where they are expected
- Phase II: Vacuum
 - Align surface to expected figure
 - ◆ Using Zernike aberrations, optimization
 - ◆ Including gravity and figure errors
 - ◆ Record all available mirror fiducial locations



Overview of Test Procedure



- Record ADM Measurement for ROC
 - Distance to null face
 - Distance to AMSD
 - Distance to Window
- Baseline Measurement Complete!
- Phase III
 - Use locations of mirror fiducials as targets for future alignments



Detail (for Be): 2 Nulls, 1 AMSD





- Be CTE results in a $\sim 13\text{mm}$ change in ROC
 - Two nulls: 10013mm and 10000mm ROC
 - ◆ 200K Measurement made with both nulls
 - ◆ Designs require 6.5mm refocus and 0.136mm horizontal realignment of pallet when switching out null
 - ◆ Views with two designs are horizontally displaced from each other (FoV's shifted).
 - + Affects fiducial targets in pixel space for final cryo measurement
 - Cryo shrinkage means fiducials have moved relative to vertex of best fit parabola.
 - ◆ Affects intermediate temperature measurements



Detail (for Be): 2 Nulls, 1 AMSD

How fiducials move



How we are aligned at
intermediate temperatures



Detail (for Be): 2 Nulls, 1 AMSD



- The transformation in parent space of ambient fiducials targets to cryo fiducial targets is required:
 - This approach means that perfect fiducial alignment at all temperatures (ROC's) for a perfect mirror would show no misalignment induced aberrations only at the initial ambient temperature (10013mm) and at the final cryo temperature (10000mm).
- Future tests should optimize differences and realignments between two null designs.



Detail:

Acquiring Mirror Fiducials

- This process is quite repeatable,
 - Charts to follow...
- but it can be time consuming, depending on number of fiducials,
- and subject to variation due to:
 - pupil focus changes
 - mirror tilt changes
 - operator changes
- So, if this approach is used in future tests
 - Never change pupil focus (not changed during AMSD tests)
 - Align to nearly same fringe pattern for every alignment
 - Use same operator OR train operators to be consistent



Detail:



Mirror Fiducial Repeatability

- Acquired 3 sets of mirror fiducials for a single alignment
- Calculated mirror alignments with all 3 sets
- Alignments varied by:
 - 0.05mm Std.Dev. Horizontal
 - 0.05mm Std.Dev. Vertical
 - 9 seconds Std.Dev. Roll
 - ♦ Equivalent to 0.11mm Vertical



Detail:



Null Fiducial Repeatability

- Acquired 4 sets of null fiducials for a single alignment
- Calculated mirror alignments with all 4 sets
- Alignments varied by:
 - 0.15mm Std.Dev. Horizontal
 - 0.09mm Std.Dev Vertical
 - 42 seconds Std.Dev Roll
 - ♦ Equivalent to 0.27mm Vertical
 - ♦ Due to instability of this roll calculation, we only use translations



Detail: Total Repeatability

- RSS of Horizontal Errors: 0.16mm
 - Corresponds to 19nm surface rms
- RSS of Vertical Errors (No Null Roll): 0.15mm
 - Corresponds to 18nm surface rms
- RSS H and V: 26nm Surface Uncertainty
- To allow the possibility of some real null roll
 - Widen RSS of Vertical Errors: 0.20mm
 - Corresponds to 22nm surface rms
- Full Surface Uncertainty: 29nm rms



Detail: IPI Repeatability

- 5 maps, each with 25 averaged frames, were acquired in sequence
- Tilt was removed in Opticode
- "Averaging" option applied in Opticode
 - Generates a report on the statistics of the sets
- Surface rms's varied by 2nm rms



Detail:

Null Installation Repeatability

- Null was removed, then reinstalled, then a Map was acquired: 5 times
 - Null alignment was not monitored / adjusted
- Tilt was removed in Opticode
- "Averaging" option applied in Opticode
- Surface rms's varied by 5nm rms
- NOTE: This test occurred over a very short time, not days.



Detail:

System Reproducibility (1)

- A good alignment was acquired.
 - 5 Maps were taken
- Entire system was significantly misaligned.
 - Null was snapped out and in 5 times
 - Hexapod was yawed back and forth 5 times
 - Hexapod was translated in Horiz & Vert
 - Mirror stage (NMSS) used to reacquire return beam
- A new good alignment was acquired and maps obtained.
 - 5 Maps were taken



Detail:

System Reproducibility (2)

- Computed misalignments of each measurement
- Alignments matched to within:
 - 0.021mm vertical
 - 0.136mm horizontal
 - 16nm rms surface
- This is within calculated uncertainties



Additional Items



- ADM measurements are extremely precise.
- Know your stages.
- Set a single coordinate system.
- Segment roll is critical to monitor and correct.
- Pupil image quality shows significant variations over the full pupil.



Summary



- Specify multi-null designs carefully.
- Improve or Remove (bypass) null marking.
- If mirror fiducial approach is used
 - NEVER change the pupil focus
 - Align to nearly same fringe pattern for every alignment
 - Use same operator or consistent operators