4D PhaseCam Capabilities


Multiple-Wavelength Mirror Phasing Modal Analysis

NASA Tech Days September, 2003

Michael North Morris, Ph.D. James Millerd, Ph.D., and John Hayes, Ph.D.



4D PhaseCam Technology



Single Frame PSI

Benefits:

- High resolution interferometric measurement
- Insensitive to vibration & turbulence
- Easy to set up and use

$$\tan \varphi = \left(\frac{B - D}{A - C}\right)$$




TWO WAVELENGTH MEASUREMENTS



Multiple Wavelength PhaseCam

• 2 Wavelengths can be used to extend range



2 Frequencies beat together to form a long equivalent wavelength

• A measurement is made at each wavelength


$$\Delta opd = 2\Delta z = \frac{\Delta \phi_e}{2\pi} \lambda_e$$

$$\Delta \phi_e = \Delta \phi_1 - \Delta \phi_2 \qquad \lambda_e = \frac{\lambda_1 \lambda_2}{|\lambda_1 - \lambda_2|}$$



Surface Subtraction

A) Surface domain (specular surfaces only)






Surface Subtraction

B) Interferogram domain (diffuse and specular surfaces)

Stetson - 8 frame phase-difference


$$\Delta Z = \frac{\lambda}{2} \arctan \left(\frac{[D_0 - B_0][A(t) - C(t)] - [A_0 - C_0][D(t) - B(t)]}{[A(t) - C(t)][A_0 - C_0] + [D_0 - B_0][D(t) - B(t)]} \right)$$






Measurement of polished and unpolished substrates

Single λ phasemap



Discontinuities not resolved
Diffuse surface - random phase
Unwrapping errors

Dual λ **phasemap**




True 3D measurement

Discontinuity resolved
1 fringe = 3 mm




2 Wavelength Measurement Uncertainty






Synthetic Wavelength Coverage




- Uncertainty at 0.115 mm must be less than 0.0014 synthetic waves
- Uncertainty at 0.018 mm must be less than 0.0088 synthetic waves
- Need three sources to cover wavelength range




Measurement Uncertainty Vs. Synthetic Wavelength for Different Source λ Uncertainties (A & C)



Measurement Uncertainty Vs. Source Uncertainty for Different Synthetic Wavelengths (A & B)



Frequency Response





Project Status

- Uncertainty Analysis is Complete
- Sources Have Been Selected
- Coupling Configuration Has Been Designed
- Optical Components Have Been Selected
- Source Module Components Are On Order
- Mechanical Design is 80% Complete
- Measurement Process Has Been Outlined
- Software Changes Are in Progress

Project is On Schedule!!



Modal Analysis



Vibrational measurement of mirrors


- •New Generation Mirrors:
 - large mirrors (10's meters)
 - lightweight
 - multi-laminate
- Susceptibility to low freq. resonant modes
 - Performance and integrity issues
 - FEM used in design
 - Actual measurements tell whole story
- Measurement requirements
 - High resolution (sub-wavelength)
 - No mass loading

Interferometry is needed



Effects of vibration

<u>Acquisition time ~ 30 microseconds</u>



Relative motion between PhaseCam and test object


Subtract two measurements to find vibration!

Phase relationship is fixed Measure surface + tilt






Synchronous Measurement




Synchronous – shutter laser

Acousto-optic modulator ~ sub-microsecond






Operational Limits






Several Resonant Modes






Phase Sweep




Al mirror, 55 Hz, first order mode




Phase Sweep







Phase Sweep



Al mirror, 3069 Hz, higher order mode



Conclusions

- PhaseCam Capabilities
 - High resolution surface figure measurement
 - Two-wavelength measurements
 - Differential surface measurement (specular and diffuse surfaces)
- Multiple Wavelength PhaseCam Development is on Schedule
- Modal Analysis Has Been Implemented
 - − Wide frequency response (sub Hz − 100kHz)
 - Accurate determination of modal deflection (<10nm)
 - Sweeping over frequency and/or phase



CONTACT INFORMATION

4D Technology Corporation

3280 E. Hemisphere Loop, Suite 112

Tucson, Arizona 85716

520-294-5600

www.4dtechnology.com

